

THE HYDRAFACIAL™ SYSTEM



The award-winning HydraFacial System advances skin health by merging invigorating spa therapies with advanced medical technology.

Treatment Boosters

Tailor your treatment with award-winning serums to address specific skin concerns



Ergonomic Handpiece

Allows control over serum type and flow

Mechanical Lymphatic Drainage

Offers further detoxification benefits

Proprietary Serums

Award winning products deliver nutrients to promote skin health



Touch Screen

Provides built-in protocols and technology updates

Skin Imaging

Enhance the client experience by demonstrating treatment results through before and after photos


Red and Blue LED Light

Allows time for product absorption and LED light benefits



System components may differ depending on the model.





Vortex-Fusion®

The unique Vortex-Fusion delivery system is what sets HydraFacial™ apart from other resurfacing procedures. HydraFacial uses a series of unique HydroPeel® Tips. Each one has a specific use to exfoliate, extract and hydrate the skin.

The spiral design of the HydroPeel Tips is used in conjunction with the HydraFacial proprietary vacuum technology and serums. This creates a vortex effect to easily dislodge and remove impurities while simultaneously introducing hydrating skin solutions with potent antioxidants.

HydraFacial Proprietary Tips



MECHANICAL LYMPHATIC THERAPY



Mechanical Lymphatic Therapy helps to detoxify the skin. Removing unwanted build-up is important in optimizing skin health so skin can be naturally healthy and radiant from the inside-out.

Treatment Areas

- Face (33% of the lymphatic system is in the face)
- Arm
- Leg
- Abdomen
- Back

*Check with your state to see who is licensed to perform mechanical lymphatic drainage.

LED LIGHT THERAPY

Light therapy is the application of light energy to the skin for therapeutic benefit. LED light therapy provides an excellent finishing touch to HydraFacial treatment.



Red LED Lights

- 12 Red LED Lights (wavelength = 620-630 nanometers)
- 12 Infrared Lights (wavelength = 800-900 nanometers)



Blue LED Lights

- Targets p. acne bacteria, the strain of bacteria that plays a key role in causing most acne conditions.